



IN THIS ISSUE

Three Tips for Faster Western Blotting

[20% off 3 Great RNA Products](#)

[Go Green with Formaldehyde-Free RNA Gels](#)

Webinar: [Go Green with Safer Laboratory Reagents](#)

Three Tips for Faster Western Blotting

Traditional Western Blotting that takes a whole day or more to complete is a thing of the past. Follow these three tips and realize significant time savings without sacrificing quality or sensitivity.

- **Tip 1.** Go from gel polymerization to visualization of your protein bands in under 45 minutes! Eliminate the stacking gel, crank up the voltage and visualize fluorescent protein bands with no post-staining and no destaining. AMRESCO's ready-to-pour [Fluorescent Sprint Next Gel™](#) acrylamide solutions give you the high performance jump start you need for time savings.
- **Tip 2.** Reduce blocking time to 10 minutes (or less) and experience enhanced chemiluminescent signals with protein-free [Rapid Block™](#). Use at least 5X less antibody to visualize the same strong signal you were used to seeing with your previous blocking buffer.
- **Tip 3.** Don't waste valuable reagents and spend valuable time running an entirely new Western to detect another protein in your sample. Strip your blot with non-hazardous, odor-free [Gentle Review™ Stripping Buffer](#). Blots can be stripped and re-probed several times without significant loss of protein.

Use Fluorescent Sprint Next Gel™, RapidBlock™ and Gentle Review™ for fast Western Blotting.



20% Discount on 3 great RNA products


[RiboZol™ RNA Extraction Reagent](#)

[RiboZol™ Plus RNA Purification Kit](#)

[Formaldehyde-Free RNA Gel Kit](#)

[More Information >](#)

Formaldehyde-Free RNA Gel



Go Green with Formaldehyde-Free RNA Gels

AMRESCO's complete kit for denaturing RNA Gels totally eliminates hazardous formaldehyde and ethidium bromide. Run RNA agarose gels *right on your bench top*, no need for a fume hood! The unique loading buffer included in the kit includes tracking dyes and a non-mutagenic fluorescent RNA dye.

[More Information >](#)

AMRESCO offers more "going green" environmentally friendly products. Register for AMRESCO's FREE "Go Green with Safer Laboratory Reagents" webinar, October 14th, 2009.

October Webinar: [Go Green with Safer Laboratory Reagents](#)

November Webinar: [Tips for Faster Western Blotting](#)

December Webinar: TBD

"I think this is a great product"
Girish Shukla, Ph.D.
Cleveland State University



[Follow us on Twitter](#)



[AMRESCO Fan Page](#)



[Join us on LinkedIn](#)



[EZ-Vision webinar on YOUTUBE](#)